

# JBJS 'watch' offers tips for reducing risk of ceramic femoral-head fractures

December 23 2015

---



A picture of fragments of a fractured ceramic femoral head. Credit: JBJS Case Connector

*JBJS Case Connector*, an online case report journal published by *The*

*Journal of Bone & Joint Surgery*, has issued a 'Watch' regarding the rare but potentially catastrophic fracture of ceramic femoral heads used in hip replacements.

Ceramic hip components are often chosen for younger patients to minimize long-term wear. Ceramic femoral head fractures arise mainly from trauma; non-compatible, damaged, or contaminated femoral head/stem taper connections; or material or manufacturing defects.

This Watch was prompted by a case [report](#) by Pomeroy *et al.* in the Dec. 23, 2015 *JBJS Case Connector*, a report by Heiner and Mahoney in the Oct. 22, 2014 *JBJS Case Connector*, and a 2014 article in *BMJ Case Reports* by Tucker and Acharya. Those three reports described patients who experienced ceramic femoral-head fractures anywhere from 10 days to 14 years after hip-replacement surgery.

Because [ceramic](#) head fractures are more likely to occur from insults during or after implantation than from manufacturing defects, the Watch includes four "golden rules" surgeons can follow to reduce the risk of these events, including making sure that the tapers on both the head and stem are compatible in all dimensions. The Watch also emphasizes the importance of patient education, during which patients should be encouraged to promptly report any and all postsurgical irregularities.

"The publication of 'Watches' helps fulfill our mission to serve the orthopaedic community," commented Marc Swiontkowski, MD, co-editor of *JBJS Case Connector*. "The 'Watch' designation may encourage the orthopaedic community to either demonstrate that these are isolated, unrelated cases or sharpen the focus further by rigorously evaluating the intervention and/or reporting related cases."

**More information:** [dx.doi.org/10.2106/JBJS.CC.15.00248](https://dx.doi.org/10.2106/JBJS.CC.15.00248)

Provided by Journal of Bone and Joint Surgery

Citation: JBJS 'watch' offers tips for reducing risk of ceramic femoral-head fractures (2015, December 23) retrieved 21 July 2023 from <https://medicalxpress.com/news/2015-12-jbjs-ceramic-femoral-head-fractures.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.