

Medtronic warns of battery problem with drug pumps

8 July 2011

(AP) -- Medical device maker Medtronic is warning patients about a rare problem with its SynchroMed drug pumps that can cause them to lose battery power and fail.

The company says it has received 55 reports in which a film formed on the pump's battery, sometimes causing a loss of power.

The pump is primarily used to treat [severe pain](#) and muscle spasticity. In some spasticity patients, the abrupt loss of pump therapy can cause serious injury and death.

[Medtronic](#) is not recalling the products and does not recommend their removal, unless they are proven to be failing.

The company said in a statement that all the reported problems involved batteries manufactured before March 2005. The company is working with U.S. regulators on a new battery design to prevent the problem.

©2011 The Associated Press. All rights reserved.
This material may not be published, broadcast, rewritten or redistributed.

APA citation: Medtronic warns of battery problem with drug pumps (2011, July 8) retrieved 7 May 2021 from <https://medicalxpress.com/news/2011-07-medtronic-battery-problem-drug.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.